




 NuaBox™

Small Rapid Nucleic Acid Detection System

Palm-size Instrument, Result in 30 minutes

Faster, Smaller, Simpler

**Sample to result takes
30 minutes or less
Positive in 10min**



**Qualitative detection in
the size of a palm
Height 8cm, width 13cm**



NuaBox™ 1A
smaller than a pad

**Just one manual step of
adding samples
Work for everyone**



**No professional operator
training required**

Testing Procedure

NuaBox™ 1A



Step 1:

Nasal Swab
provided in the kit



Step 2:

Rotate Swab
dipped into tube



Step 3:

Add sample
lyophilized power inside



Step 4:

Initiate Detection
inserted into NuaBox

Two Models Available

One sample or ten

**One test at a time,
10 to 20 tests per day per machine**



 NuaBox™ 1A

**Ten tests at a time(one test, one tube)
100 to 200 tests per day per machine**




 NuaBox™ 10B

Two Single-gene-detected Reagents

ORF1ab
Higher specificity



N Gene
Higher sensitivity




- SARS-CoV-2 ORF1ab gene detection (20 tests/box [PGA4132P4])
- SARS-CoV-2 N gene detection (20 tests/box [PGA4133P4])

No Contamination

Special Patented Tube Design

- *Double Rubber Ring to Prevent Pollution*



- *Hook on the Cap to Lock the Tube Tightly*

- *Negative Pressure Cabin to Contain the Escaping Aerosol*

Much More Reliable
than Traditional Tubes



Exogenous Internal Control

NuaBox™ contains an internal control that has been designed to control for sample inhibition, amplification, and assay reagent function.

Exogenous internal control provided in the NuaBox™ lyophilized kit can exclude the invalid result among test.


Direct Isothermal Fluorescence Amplification

1. Isothermal Real-time Fluorescence Probing is a method both fast as LAMP and specific as fluorescence probing.
2. Lyophilization adopted is convenient for transportation and room temperature storage for non-lab environment.
3. Built-in professional optical parts for ensuring the sensitivity during the isothermal amplification.
4. Sensitivity: < 500 copies/mL



Nucleic Acid Amplification Test Using Anterior Nasal Swabs

Realtime Fluorescence Signal



Start testing



Temperature
59.24/61



History
8

Name P1



Sample Type Oral swab

Report

Test Item SARS2


Sample ID N4

Result Positive

Generate Report

TE GEN 透景

10000 copies/mL



Start testing



Temperature
59.07/61



History
9

Name P3



Sample Type Oral swab

Report

Test Item SARS2

Sample ID N3

Result Positive


Generate Report

TE GEN 透景

1000 copies/mL

Nucleic Acid Amplification Test Using Anterior Nasal Swabs

Realtime Fluorescence Signal



Start testing

Temperature 59.27/61

History 11

Name P6

Sample Type Oral swab

Test Item SARS2


Sample ID N2

Result Positive

Generate Report

TE GEN 透景

500 copies/mL



Start testing

Temperature 59.26/61

History 4

Name P7

Sample Type Oral swab

Test Item SARS2

Sample ID NC

Result Negative


Generate Report

TE GEN 透景


Negative

NuaBox™ Comparison Between Tellgen DIFA and Tradition LAMP

System	Detection	Observation	Sensitivity	Specificity	Pollution	Report
Tellgen NuaBox System	Fluorescence Target Probe	<i>Realtime Fluorescence Signal</i>	High	High	No Contamination (Level 2)	Digital Report
Traditional LAMP	Hydroxy Naphthol Blue	Naked Eye Discoloration	Medium	Low	5% Leakage After Amplification	Handwritten Report



Tellgen Direct Isothermal Fluorescence Amplification(DIFA)



Traditional LAMP



CE Certified Tellgen NuaBox™ Small Rapid Nucleic Acid Detection System

NuaBox™ Applicable Scenarios



E&A Dpmt

Fast result with qualitative report in minimum 10 minutes.



Emergency

Before boarding or before exhibiting in need for test result



Corporates

Places where conditions are restricted, like factories and construction sites



Others

Clinics, third party labs, schools, and nursing houses in global markets



Tellgen NuaBox™

Rapid Nucleic Acid Test, Anytime Anywhere