



## SARS-CoV-2 Nucleic Acid Detection kit Based on RealTime PCR Platform

### Stable

At room temperature for 9 months

### Cost effective

No requirements for maintaining cold temperature, save freight of both domestic and international transportation

### Economic

Save freight of dry ice shipments in air transport

## Advantages

### High specificity


A dual-target assay for the conserved region of spike protein S gene (FAM) and nucleocapsid protein N gene (ROX). Human  $\beta$ -actin gene (Cy5) acts as an internal reference control the PCR reaction.

### Stable

Lyophilized form of the reagent can be transport at room temperature.

### Multiple Applicable instrument

Any real-time PCR instrument with FAM/ ROX/ Cy5 detection channel.



### Tellgen Corporation.

Address: Building 1, No.115, Lane 572, Bibo Road, Free Trade Zone, Pudong District, Shanghai, China  
Mail: [igenhadeer@tellgen.com](mailto:igenhadeer@tellgen.com)  
Tel: +86 13600425109